



POLISH CENTRE FOR TESTING AND CERTIFICATION

02-699 Warszawa, ul. Kłobucka 23A

Mechanical Laboratory

tel.: (+48 22) 46 45 594; fax: (+48 22) 46 45 563

e-mail: labmech@pcbc.gov.pl

Date: 21.01.2010

TEST REPORT NO. BR – 241/L-496/2009

Page 1 of 4

Subject of testing: *Foldable shower chair*

Classification according to PN-EN ISO 9999:2001 :
09 33 03

Manufacturer : *MOBILEX A/S, Noerskovvej I
DK - 8660 Skanderborg*

Catalogue no.: *9041190205*

Number of specimens : *1*

Applicant: *A-Net s.c. 93-469 Łódź,
ul. Łaskowice174*

Kind of testing *Testing scope according to application of Client
Mechanical testing for conformity with PN-EN 12182:2005
excluding clause 4.1, 4.3, 5.2, 22, 24*

Test started: *10.12.2009*

Test finished : *21.01.2010*

Performed by:

Mirosław Szymański

Checked by :

Ireneusz Czerwiński

Approved by :

KIEROWNIK
LABORATORIUM MECHANICZNEGO

mgr inż. Andrzej Tkaczyk

Special comments / enclosures:

Copyright © 2005 by Polish Centre for Testing and Certification (applicable to report form)

Test results refer only to tested units.

Test results reported here are not applicable to the further modifications of the product affecting its structure, material or technology.

This test report shall be neither copied differently as in the whole nor be published without written consent of the Laboratory.





TESTING

NORMATIVE REFERENCES

PN-EN 12182:2005

Technical aids for disabled persons – General requirements and test methods

PHOTO OF PRODUCT



CHARAKTERISTICS OF PRODUCT

Name of product: <i>Foldable shower chair</i>	Maximum weight of user: 100 kg	Mass of product: 4,05 kg
--	---------------------------------------	---------------------------------

Dimensions:

Overall height:	800 mm – 900 mm
Maximum width:	555 mm
Maximum length:	510 mm
Seat height:	465 mm – 565 mm (4 levels of adjustment)
Shipping measurements:	800 mm x 555 mm x 160 mm

Material:
Aluminum alloy, plastic

