



CENTRE FOR TESTING AND CERTIFICATION - MECH-TEST

Mechanical Laboratory

05-077 Warszawa-Wesoła, ul. Klonowa 22
tel.: +48 603 23-26-45, e-mail: cbc.mech.test@gmail.com, www.cbc.org.pl

Date 22.10.2014


TEST REPORT NO. **CBC-082/2014**

Page 1 of 8

| | | |
|----------------------------|---|---|
| Subject of testing: | <i>Shower chair</i> | Classification according to PN-EN ISO 9999:2011: 09 33 03 |
| Type / Model: | <i>KAKADU</i> | SN.: 13127950087 Art. no: 302018 |
| Manufacturer: | <i>MOBILEX A/S, Grønlandsvej 5, DK-8660 Skanderborg</i> | Number of specimens: 1 |
| Applicant: | <i>A-Net s.c. 93-469 Łódź, ul. Łaskowice 174</i> | |
| Kind of testing | <i>Testing scope according to application of Client Mechanical testing for conformity with PN-EN 12182:2012, PN-EN 12183:2010 cl.7.3, 7.4.1, 7.5, 7.8, 7.10, ISO 7176-1:1999, ISO 7176-3:2003, PN-ISO 7176-8:2002 cl. 8.4, 8.5, 10.4.2, PN-EN 1021:2007</i> | |
| Test started: | 6.10.2014 | |
| Test finished: | 22.10.2014 | |

Approved by:

DYREKTOR


mgr inż. Andrzej Tkaczyk

Special comments / enclosures:

- 1) Annex 1 – Identification of wheelchair elements

Copyright © 2012 by Centre for Testing and Certification (applicable to report form)

Test results refer only to tested units.

Test results reported here are not applicable to the further modifications of the product affecting its structure, material or technology.

This test report shall be neither copied differently as in the whole nor be published without written consent of the Laboratory.



CHARACTERISTIC OF PRODUCT

| Name of toilet chair: <i>Shower chair KAKADU</i> | | SN 13127950087 Art.no 302018 | |
|--|---|---|--|
| Maximum load capacity: <i>150 kg</i> | | Overall mass of wheelchair: <i>12,75 kg</i> | |
| Description | | | Comments |
| Dimensions: | Length (footboard in maximum extension position): | | 915mm |
| | Height : | | 967mm |
| | Width: | | 547mm |
| Construction of frame: | Material: | | Aluminum alloy |
| | Method of fastening frame elements: | | Welding |
| | Folding/unfolding: | | Unfolding |
| Castor wheels | Ø of wheel: | | 125mm |
| | Width: | | 31mm |
| | Material of ring of a wheel: | | Plastic |
| | Material of fork: | | Plastic |
| | Vertical adjustment (number of fixing positions) | | NO |
| | Horizontal adjustment (number of fixing positions): | | NO |
| | Adjustment of axis inclination angle: | | NO |
| Backrest | Folding/unfolding: | | Unfolding |
| | Backrest inclination adjustment | stepless: | NO |
| | | number of fixing positions | NO |
| Tilt levers | Two singular: | | NO |
| | One lateral: | | NO |
| Push handles | Kind: | | One transverse handle |
| Parking brake | Left: | | YES |
| | Right: | | YES |
| | Kind: | | Lever , foot-push brake, wheel-mounted |
| | Material of lever: | | Plastic |
| | Fastening to frame: | | Non-adjustable brake |
| Upholstery | Material: | | Soft Plastic |
| | Colour: | | Black |
| Wheel space in forward direction position: | | | 467mm |
| Wheel space in backward direction position: | | | 467mm |
| Legrests | Common for both legs: | | NO |
| | Separate for each leg: | | YES |
| | Stationary: | | NO |
| | Folding: | | YES |
| | Vertical adjustment (number of fixing positions) | | YES |
| | Horizontal adjustment (number of fixing positions): | | NO |
| | Angle adjustment (number of fixing positions): | | NO |
| | Material of legrest: | | Aluminum alloy Plastic |
| Accessories | Seat belt | | NO |
| | Service : | | YES |

PHOTO OF PRODUCT



| | | |
|--|-------------------------------|---|
|  | Mobilex | Shower chair |
| MOBILEX A/S Grønlandsvej 5 DK-8660 Skanderborg Tel: +45 87 93 22 20 www.mobilex.dk | |  ↓ Max. 150 kg |
|  | Produced <input type="text"/> | Serial no. <input type="text"/> |