

## CENTRE FOR TESTING AND CERTIFICATION - MECH-TEST **Mechanical Laboratory**

05-077 Warszawa-Wesoła, ul. Klonowa 22 tel.: +48 603 23-26-45, e-mail: cbc.mech.test@gmail.com, www.cbc.org.pl

Date 18.05.2017

## TEST REPORT NO. CBC -086/2017

Page 1 of 4

Subject of testing:

Walking aids with built-in handgrips and three or

Classification according to PN-EN ISO 9999:2011:12 06 06

more legs of which two or more are having wheels.

which provide support whilst walking

Type / Model:

Carbon Rollator

Athlon 550 SL

Nr kat.: CRBKL700

Number of specimens: 1

SN: (01)05901912635449 (11)170601(21)0001

Manufacturer:

REHASENSE Sp. z o.o.

Sulejowska 45

97-300 Piotrków Trybunalski

Applicant:

A-Net s.c.

93-469 Łódź,

ul. Łaskowice174

Kind of testing

Mechanical testing for conformity with PN-EN ISO 11199-2: 2005

excluding clause 4.10, 6.2, 6.3

Test started: 9.05.2017

Test finished: 18.05.2017

Approved by:

Special comments / enclosures:

Copyright © 2012 by Centre for Testing and Certification (applicable to report form)

Test results refer only to tested units.

Test results reported here are not applicable to the further modifications of the product affecting its structure, material or technology.

This test report shall be neither copied differently as in the whole nor be published without written consent of the Laboratory.



Page: 2 of 4

## CHARACTERISTIC OF PRODUCT

Name: Carbon Rollator Athlon 550 SL

Dimension of rollator:

Nr. kat: CRBKL700

SN: (01)05901912635449(11)170601(21)0001

Maximum permissible user mass: 200 kg

Mass of rollator: 6,265 kg

РНОТО	OF	PRO	OD	UCT

11.252.2	Descript		5
E	Elements/parameters/materials/		Comments
	Distance between	560 mm	
L L	handgrips (dimension 2)		
Dimensions od walking rollator (fig. 2 PN-EN ISO 11199-2)	Angle between of handgrip	00	
rol 99-	axis and direction of		
ng 	movement (α)		
is o	Height of rollator	743/770 mm	min.
ws VIS	(dimension 6)	1015 mm	max.
P F	Width of rollator	682 mm	
Suc PN	(dimension 5)		
ısic		Turning width 890 mm	
ner (fig	(dimension 1)		
D.	Length of rollator 670mm		
	(dimension 4)		
Dimen	sions of folded rollator (mm)	855 x 660 x 22	8
33	Handgrip - diameter	37 mm	Anatomical
à		1.0	handgrip
	Handgrip - length	120 mm	
	Front wheels- quantity	2	castor
Wheels of rollator	Front wheels - diameter	202,5 mm	wheels
llo.	Front wheels – width		
<u>ب</u>	Front wheels - brake	none	
ls o	Rear wheels - quantity	2	
iee	Rear wheels - diameter	202,5 mm	]
₩	Rear wheels - width	34 mm	
**	Rear wheels - brake	Included	
Tip	Diameter	NAME OF A STATE OF THE STATE OF	DO. VANDO DO DO
	Material	Not any	
	Colour		
4	Front legs	Carbon,	
fro	Bracing member (no. 8)	Hard plastic,	
	Rear legs	Bolts, nuts	
al c	Height adjusting device (no.		
eri	4)		
Material of tor (fig. 1	Handgrip (no 5),	Hard plastic	
Z +3	Brake elements		

















Report no.: CBC-086/2017

Page: 3 of 4

		RES	ULT	OF TESTS ACCORDIN	NG TO PN-E	V ISO 1110	0 2.20	1 age. 3 0j 4
		TES	OLI	Of TESTS ACCORDIT	NG TO TIN-E	130 1119	T	05
Requirements according to clause	Test me- thod according to			Checked naracteristics/assemblies/parameters		Real	Test result	Comments
4.1	Measur.		Manoeuvrability			ø 202,5 mm width 34mm Conf.	Pos.	ø front wheels ≥75mm outdoor intended rollator ø front wheels ≥180mm width of wheels ≥28mm
4.2	5.3	Forw	ard-dire	ection stability		17,0° Conf.	Pos.	Stability required ≥ 10°
4.2	5.4 5.5			rection stability		8,5 ° Conf.	Pos.	Stability required ≥ 7°
4.2	5.6	Stabil		ection stability	forwards	4,0 ° Conf. 15,8 ° Conf.	Pos. Pos.	Stability required $\geq 3.5^{\circ}$ Stability required $\geq 10^{\circ}$
				basket, bag, drip, oxygen cylinder	backwards	13,3 ° Conf.	Pos.	Stability required ≥ 7°  Stability required ≥ 7°
	S		side		side	6,0° Conf.	Pos.	Stability required ≥ 3,5°
4.3	V/I	V	vheels	g facility during rollator motion w	N V 1889	Conf.	Pos.	
	V/I	r	esting s	brakes in rollator with more than a eat or intended for outdoor use	2 wheels and	Conf.	Pos.	
	5.7.1.1			rip distance (fig. 4, dimension 1)		65 mm Conf.	Pos.	≤75 mm
	5.7.1			brake effectiveness		Conf.	Pos.	Movement of rollator ≤10 mm in 1 minute
	Measur.			set parking brake		20N Conf.	Pos.	≤ 60 N
	Measur.	I		release parking brake		10N Conf.	Pos.	≤ 40 N
	5.7.2			brake effectiveness		Conf.	Pos.	Movement of rollator ≤ 10 mm in 1 minute
	V/I			ty to compensate brake wear		Conf.	Pos.	
	V/I			ot adversely affected by folding, ung actions of rollator	nforlding or	Conf.	Pos.	
4.4	Measur. V/I	Handg	grip			37 mm Conf.	Pos.	Width of handgrip ≥20mm and ≤50mm
4.5	Measur. V/I	Leg se	ection a	nd tip			N/A	ø tip ≥35mm (tested rollato is equipped in four wheels)
4.6	5.10		T3	static loading durability		Conf.	Pos.	1 minute under load 1,2 x user's weight±2% (240kg)
4.7	5.12		hanical ability			Conf.	Pos.	200 000 cycles with load. 160kg±2%, f=1Hz
4.7	5.11			Static loading test		Conf.	Pos.	loading 240kg±2%, 5sek. NOTE 1
4.8	V/I		ing dev			Conf.	Pos.	
4.9	5.14	Foldin	g mech			Conf.	Pos.	
4.11	ISO 10993-	uls sh	В	iocompatibility of material with h	uman body		N/T	
	V/I	Materials and finish		ree of discolouring of skin or cloth rith rollator materials	ning in contact	Conf.	Pos.	
	V/I	≥ ≅	В	urrs, shar edges, projections		Conf.	Pos.	
<u> </u>			,	Marking and labe	lling of product			
6.2	V/I	a) Max	imum ı	iser mass		Conf.	Pos.	
	H	c) Max	imum s	safe working load (SWL) to be ma allowed angle between the longitude	rked on accessories	Conf.	Pos.	
			and the	e direction of motion, if the handle		Conf.	Pos.	
		d) Man	) Manufacturer's name or trade name and address			Conf.	Pos.	
				er's model identification name and	d/or number	Conf.	Pos.	
				year of manufacture		Conf.	Pos.	
			imum e	extension of the height adjustment.	, marked on the	Conf.	Pos.	
		h) Max	imum v	width of the rollator		Conf.	Pos.	
		78		ended for outdoor/indoor use		Conf.	Pos.	
1.10			Varning showing allowed angle between handle axis and direction ovement or physical stop of angle adjusting			of	N/A	



Mechanical Laboratory of CBC		Report	Report no.: CBC-086/2017			
			Î	I	Page: 4 of 4	
		Contents of user manual and/or assembly manual or clear and indeli	ole markin	g of produ	et	
6.3	V/I	a) Maximum rollator height	Conf.	Pos.		
		b) Minimum rollator height	Conf.	Pos.		
		c) maintenance and cleaning instructions, including a description of the method and suitable cleaning agents and any precautions needed to avoid corrosion and/or ageing of the materials used in construction of the rollator	Conf.	Pos.		
		d) Instructions for assembly, adjustment of all kinds, folding and unfolding	Conf.	Pos.		
		<ul> <li>e) Warnings and advice about precautions relating to safe distances between moving and stationary parts (see EN 12182, Clauses 12 and 13, for guidance)</li> </ul>	Conf.	Pos.		
		f) Maximum safe working load (SWL) for load carrying accessories such as basket, tray, shopping bag, etc.	Conf.	Pos.		
4.10	V/I	Warning in user manual on consequences of such an adjustment of angle between handle longitudinal axis and direction of movement outside allowed value (when handles are adjustable aside).		N/A		
		TEST CONDITIONS				
Ambient temperature			19°C		Required temperature 21°C ±5°C	
Relative humidity of air:			55 %		Not required	
Commen	701					
		n maximum height adjustment of rollator.				
All tests p	performend in	the least stabble position of self-adjusting wheels.				
ests per	formed with ha	andles positioned at their maximum (allowed) angle to the direction of motion	on (when ac	ljustment is	s possible).	
equence	of tests: stabil	lity test, static loading test, fatigue test.				
	tor was tested.					
During v	isual inspect	ion before testing any visible defects that could have influence on te	st results v	vere not st	ated.	
		g = pegative: N/T = pot tested: N/A = pot applicable: N/B = pot require				

Pos. – positive; Neg – negative; N/T – not tested; N/A – not applicable; N/R – not required , N/O – not occurred , V/I.- visual inspection, Conf.- conformed.

NOTE 1: Deformation - 68mm, elastic deformation - 64 mm, permanent deformation - 4mm (0,4%).

## **CONCLUSIONS:**

Testing object conforms with requirements of PN-EN ISO 11199-2: 2005, excluding clause 4.10, 6.2, 6.3 in scope of mechanical testing ordered by client, excluding biocompatibility tests of material with human body according to PN-EN ISO 10993-1:2010.

---- END ----

